



**Azadeh Hadadi**

**Report Concepts of Medical Image  
Post-Processing (Management and  
Post-Processing of Prostate MRI)**

**15-May-2020**

## 2020 MSc. COMPUTER VISION

I

Reporter (name/surname):

**Azadeh Hadadi**

Module (category):

**Medical imaging**

Report Title:

**Report Concepts of Medical Image Post-Processing**

Sub-title:

**Management and Post-Processing of Prostate MRI  
(Programing in MATLAB)**

Professor/lecturer/instructor:


**Professor Christian Mata Miquel**

Date:

**15-May-2020**

At:

**Centre université Condorcet, Le Creusot, France**



## TABLE OF CONTENTS

TABLE OF FIGURES .....	4
1. ABSTRACT.....	6
2. INTRODUCTION.....	6
2.1 Problem Definition.....	6
2.2 Background .....	6
2.3 Report Organization .....	6
3. Prostate structure .....	7
4. Prostate Cancer.....	7
5. Prostate cancer detection using MRI .....	8
5.2 ANATOMICAL imaging– 3D-T2 weighted imaging .....	9
6. principle of Snake Contour .....	9
6.1 SNAKE Behavior .....	9
6.2 Image fources .....	10
6.3 Line Functional.....	10
6.4 Edge Functional .....	10
6.5 Scale Space .....	10
6.6 Termination Functional.....	10
6.7 Constraint Energy .....	11
6.8 Optimization through gradient descent .....	11
6.9 Discrete approximation .....	11
7. Procedure.....	12
7.1 Actions .....	12
7.2 Snake Contour .....	13
7.3 3D Representation .....	15
8. Flow charts and algorithms.....	15
8.1 Function flow charts .....	15
9. Results and functional display .....	20

10.	CONCLUSION .....	23
11.	REFERENCES .....	24
12.	Appendix.....	24
12.1	loadDICOM .....	24
12.2	Anonymize .....	24
12.3	DICOMToJPG.....	24
12.4	JPGToDICOM.....	25
12.5	Smooth .....	25
12.6	Snake .....	25
12.7	CalculateSurface .....	25
12.8	CalculateVolume.....	25
12.9	3D Representation .....	25

## TABLE OF FIGURES

Figure 1: different zones in a prostate .....	7
Figure 2:50-year-old male with prostate cancer. (a) Axial T2WI showing low signal within the rightward aspect of the peripheral zone.Note the normal-appearing contralateral peripheral zone (PZ) comprised of glandular elements. (b) Axial T1WI showing an enlarged rightexternal iliac chain lymph node.....	8
Figure 3: actions panel: Load DICOM and JPG to DICOM buttons are active .....	12
Figure 4: successful anonymized action message .....	12
Figure 5: successful DICOM to JPG message .....	13
Figure 6: error message in converting JPG to DICOM .....	13
Figure 7: successful JPG to DICOM message .....	13
Figure 8: sigma value for Gaussian filter.....	13
Figure 9: snake panel.....	14
Figure 10: snake contour .....	15
Figure 11: 3D representation panel.....	15
Figure 12: warning folder out of JPG file .....	15
Figure 13: load DICOM function's flow .....	16
Figure 14: Anonymize function's flow .....	16

Figure 15: DICOM to JPG function's flow .....	17
Figure 16: JPG to DICOM function's flow .....	17
Figure 17: smooth function's flow .....	18
Figure 18: Snake function's flow .....	18
Figure 19: (a) calculate Surface's function flow (b) calculate volume's function flow .....	19
Figure 20: 3D representation function 's flow .....	19
Figure 21: DICOM files .....	20
Figure 22: status panel .....	20
Figure 23: anonymized DICOM file .....	20
Figure 24: JPG and mat files .....	21
Figure 25: smoothed image .....	21
Figure 26: segmented DICOM file (a) peripheral zone (b) transitional zone (c) tumour (d) central zone .....	21
Figure 27: snake parameters .....	22
Figure 28: surface and volume of each zone .....	22
Figure 29: created masks .....	22
Figure 30: all zones in one image .....	23
Figure 31: 3D representation of zones .....	23
Figure 32: 3D representation of full prostate gland .....	23

## 1. ABSTRACT

This report will present a management and post-processing tool for prostate MRI analysis. The objective of the tool and requested functionalities will be presented and usage procedure will be detailed. The application has a user interface (GUI) developed in MATLAB code and is easy to use. This GUI application is all about different operations related to DICOM format manipulation such as extraction patient's information, anonymization, converting to JPG files and vice versa, prostate zone segmentation through Active Contour Model (Snakes), calculating the surface and volume of each zone, and finally, show segmentation result in a 3D representation. A Database including DICOM image of 64 patients is the input data to the tool and the result will be displayed in the GUI interface. At the end of the report some images will be processed and the result will be presented. The report will end up to the conclusion and a summary of the findings of this work.

## 2. INTRODUCTION

### 2.1 PROBLEM DEFINITION

The objective of the development is to design and implement a GUI interface in MATLAB capable of performing following functionalities on a set of prostate images with DICOM format:

#### "A. 1st stage

- Display the information of DICOM (PatientName, PatientID, PatientBirthDate, StudyID, StudyDate, SliceLocation, and InstanceNumber)
- Anonymize DICOM images: Modify the following fields (PatientName, PatientID, and BirthDate) for all the set of images
- Save the anonymized images (Create a new folder containing the new anonymized images)
- Convert the DICOM image in JPG format (Be careful with the DICOM information!)
- Convert the JPG image into DICOM format

#### B. 2nd stage

- Snake segmentation for each region (only ZP, ZT, ZC, and the tumour region)
- Show a 3D representation of the prostate gland. (All the prostates)
- Show a 3D representation for PZ and CZ
- Show a 3D representation of the tumour region.
- \* For the 3D representation should do using a transparent system.

#### C. 3rd stage

- for each region calculate the surface (Surface = number of pixels x spatial resolution)
- Calculate the volume for each region."

### 2.2 BACKGROUND

The prostate is a small walnut shaped gland in the pelvis of men. Prostate cancer is the development of cancer in the prostate, a gland in the male reproductive system. In the last few decades, new imaging techniques based on Magnetic Resonance Imaging (MRI) have been developed to improve diagnosis. In practice, diagnosis can be affected by multiple factors such as observer variability and visibility and complexity of the lesions. In this regard, computer-aided detection and computer-aided diagnosis systems have been designed to help radiologists in their clinical practice. Research on computer-aided systems specifically focused for prostate cancer is a young technology and has been part of a dynamic field of research during the last 10 years.

### 2.3 REPORT ORGANIZATION

The remaining of the report is organized as follows. In section 3, the prostate structure and different zones will be detailed. Section 4 will discuss prostate cancer. In section 5, prostate MRI images and their detail will be explained. Thereafter, the theoretical aspects of the snake algorithm proposed by Kass et al. [1] will be explained and theoretical background and mathematical equations will be introduced in more detail in section



6. The application procedure will be detailed in section 7. This section will briefly explain how a new user shall install and use the GUI to process a given patient prostate MRI images. All of the steps and as per above functionality (see section 2.1 problem definition sub-section) will be described in this section. Section 8 explains all algorithms used in the development with the main functions' flow charts. Section 9 contains the results obtained running the application with comparison with some other similar literatures. The report will end up with conclusion and references. Some parts of the main functions is attached in the Appendix I. the complete functions can be found in the attached MATLAB code (see ZIP file attached).

### 3. PROSTATE STRUCTURE

The prostate is a walnut-sized gland located between the bladder and the reproduction organs of the males. The prostate is situated right in front of the rectum. The urethra runs through the center of the prostate, from the bladder to the penis, letting urine flow out of the body. The prostate secretes fluid that nourishes and protects sperm. During ejaculation, the prostate squeezes this fluid into the urethra, and it's expelled with sperm as semen. The vasa deferentia (singular: vas deferens) bring sperm from the testes to the seminal vesicles. The seminal vesicles contribute fluid to semen during ejaculation.

In general, doctors divide the prostate into three zones (see Figure 1below): the peripheral (1), central (2), and transition (3) zones. Most prostate cancers arise in the peripheral zone, the outer area of the prostate, which is next to the rectum. Doctors can sometimes feel a prostate tumor through the rectal wall during a digital rectal examination. Semen stored in the seminal vesicles passes through the central zone and into the urethra for ejaculation. The transition zone is the innermost section of the prostate. Prostate tissue in this area wraps around the urethra.




Figure 1: different zones in a prostate

### 4. PROSTATE CANCER

Prostate cancer is the second most diagnosed cancer of men all over the world. In France, prostate cancer is the most common male cancer (71,200 new cases estimated in 2011) and the third leading cause of cancer death in men (8,700 deaths per year). The median age of diagnosis is 74 years, and 44 Percentage of prostate cancers are diagnosed after 75 years. The average age of death for prostate cancer is 78, almost the average life expectancy of men in France.

In the last few decades, new imaging techniques based on Magnetic Resonance Imaging (MRI) have been developed to improve diagnosis. In practice, diagnosis can be affected by multiple factors such as observer variability and visibility and complexity of the lesions. In this regard, computer-aided detection and computer-aided diagnosis systems have been designed to help radiologists in their clinical practice. Research on computer-aided systems specifically focused on prostate cancer is a young technology and has been part of a dynamic field of research for the last 10 years.

Prostate cancer typically presents as a round or ill-defined low signal-intensity focus in the peripheral zone on T2WI. Since the majority of prostate carcinomas arise in the peripheral zone, many of them can be readily detected within the high signal-intensity background of the loosely packed normal peripheral zone glandular tissue as is shown in Figure 2[1].




Figure 2: 50-year-old male with prostate cancer. (a) Axial T2WI showing low signal within the rightward aspect of the peripheral zone. Note the normal-appearing contralateral peripheral zone (PZ) comprised of glandular elements. (b) Axial T1WI showing an enlarged right external iliac chain lymph node.

## 5. PROSTATE CANCER DETECTION USING MRI

Magnetic resonance imaging (MRI) is a medical imaging technique used in radiology to form pictures of the anatomy and the physiological processes of the body which plays an important role in the anatomic evaluation, detection, and staging of prostate cancer is well established.

Multiparametric MRI (mpMRI) represents a growing modality for the non-invasive evaluation of prostate cancer, PCA and is increasingly being used for patients with persistently elevated PSA and prior negative biopsies, for monitoring patients in active surveillance protocols, for preoperative characterization of cancer for surgical planning, and in planning for MRI-targeted biopsy.

Prostate MRI is performed using either 1.5-T or 3-T magnetic field strengths, typically with the combined use of endorectal and pelvic phased-array coils to maximize the signal-to-noise ratio. A bowel relaxant will also optimize the study by reducing artifact from bowel motion. The mpMRI is the current reference standard because no single MRI sequence is entirely sufficient to characterize prostate cancer. The optimal combination and interpretation approach of anatomic and functional MR sequences still needs to be established. However, the more functional sequences that are combined, the better the accuracy appears to be. Recently, Turkbey et al. [ ] reported that a four-sequence multiparametric MRI (T2-weighted imaging, DWI, DCE-MRI, and MRS) had sensitivity of 86% and specificity of nearly 100% in a prospective trial of 45 patients. A number of studies that evaluated the use of a four-sequence multiparametric MRI approach in the diagnosis of localized prostate cancer reported sensitivity, specificity, accuracy, PPV, and NPV for the detection of prostate cancer of 69–95%, 63–96%, 68–92%, 75–86%, and 80–95%, respectively.

The European Society of Urogenital Radiology (ESUR) and the European Association of Urology (EAU) have recently published clinical guidelines for multiparametric MRI of prostate outlining both minimal and optimal requirements to allow a more consistent and standardized approach. Both articles recommend including T1-weighted, T2-weighted, DWI, and DCE-MRI sequences, but the addition of MRS is optional. The ESUR guidelines also outline the prostate imaging reporting and data system (PI-RADS) structured reporting system, which includes a 5-point scale for reporting the likelihood of clinically significant prostate cancer and probability of extraprostatic disease being present. The value of PI-RADS as a diagnostic tool and as a predictor of patient outcomes remains to be determined [4].



## 5.2 ANATOMICAL IMAGING– 3D-T2 WEIGHTED IMAGING

T2-weighted imaging High-resolution axial, sagittal, and coronal T2-weighted imaging (T2WI) sequences offer excellent soft tissue contrast and depiction of the zonal anatomy of the prostate. As such T2WI is best placed to identify defects in zonal anatomy, as exhibited by PCA cells (orto depict seminal vesicle invasion and extra-capsular extension of disease). Normal PZ tissue is extremely water rich and composed of numerous ductal and acinar elements with sparsely interwoven smooth muscle. This gives it a bright or high-intensity appearance on T2-weighted images. On the other hand, PCA in the PZ appears as a rounded or ill-defined low-signal intensity focus that contrasts with the high-signal intensity of the loosely packed normal PZ tissue.<sup>51</sup> The presence of a low-signal intensity focus in the PZ of T2-weighted images does not definitively indicate the presence of cancer because conditions such as prostatitis, atrophy, and prior biopsy-related hemorrhages may mimic this appearance as well. The normal TZ tissue has less water content, more compact smooth muscle, and sparser glandular components than the PZ and thus appears relatively darker on T2WI. PCA in the TZ appears as a homogeneous mass of low-signal intensity with indistinct margins. It is often hard to distinguish cancer from stromal BPH in the TZ, which may also appear as low-signal intensity due to its high muscular and fibrous contents [4].

## 6. PRINCIPLE OF SNAKE CONTOUR

A snake is an energy-minimizing spline guided by external constraint forces and influenced by image forces that pull it toward features such as lines and edges. Snakes are active contour models: they lock onto nearby edges, localizing them accurately. Scale-space continuation can be used to enlarge the capture region surrounding a feature. Snakes provide a unified account of a number of visual problems, including detection of edges, lines, and subjective contours; motion tracking; and stereo matching [2].

### 6.1 SNAKE BEHAVIOR

Basic snake model is a controlled continuity spline under the influence of image forces and external constraint forces. The internal spline forces serve to impose a piecewise smoothness constraint. The image forces push the snake toward salient image features like lines, edges, and subjective contours. The external constraint forces are responsible for putting the snake near the desired local minimum. These forces can, for example, come from a user interface, automatic attention mechanisms, or high-level interpretations.

Representing the position of a snake parametrically by  $v(s) = (x(s), y(s))$ , we can write its energy functional as:

$$E_{snake}^* = \int_0^1 E_{snake}(v(s)) ds = \int_0^1 (E_{internal}(v(s)) + E_{image}(v(s)) + E_{con}(v(s))) ds$$

Where  $E_{int}$  represent the internal energy of the spline due to bending,  $E_{image}$  gives rise to the image forces, and  $E_{con}$  gives rise to the external constraint forces. The internal spline energy can be written

$$E_{internal} = E_{cont} + E_{curv}$$

$$E_{internal} = \frac{1}{2}(\alpha(s) |v_s(s)|^2) + \frac{1}{2}(\beta(s) |v_{ss}(s)|^2) = \frac{1}{2} \left( \alpha(s) \left\| \frac{d\bar{v}}{ds}(s) \right\|^2 + \beta(s) \left\| \frac{d^2\bar{v}}{ds^2}(s) \right\|^2 \right)$$

The spline energy is composed of a first-order term controlled by  $\alpha(s)$  and a second-order term controlled by  $\beta(s)$ . The first-order term makes the snake act like a membrane and the second-order term makes it act like a thin plate. Adjusting the weights  $\alpha(s)$  and  $\beta(s)$  controls the relative importance of the membrane and thin-plate terms. Setting  $\beta(s)$  to zero at a point allows the snake to become second-order discontinuous and develop a corner [2].

## 6.2 IMAGE FOURCES

In order to make snakes useful for early vision we need an energy functional that attracts them to salient features in images. The total image energy can be expressed as a weighted combination of the three energy functions:

$$E_{image} = w_{line}E_{line} + w_{edge}E_{edge} + w_{term}E_{term}$$

By adjusting the weights, a wide range of snake behaviour can be created [2].

## 6.3 LINE FUNCTIONAL

The simplest useful image functional is the image intensity itself. If we set

$$E_{line} = I(x, y)$$

Then, depending on the sign of  $w_{line}$ , the snake will be attracted either to light lines or dark lines. Subject to its other constraints, the snake will try to align itself with the lightest or darkest nearby contour [2].

## 6.4 EDGE FUNCTIONAL

Finding edges in an image can also be done with a very simple energy functional. If we set

$$E_{edge} = -|\nabla(I(x, y)^2)|$$

Then the snake is attracted to contours with large image gradients.

## 6.5 SCALE SPACE

A snake originating far from the desired object contour may erroneously converge to some local minimum. Scale space continuation can be used in order to avoid these local minima. This can be achieved by a blurring filter on the image and reducing the amount of blurring as the calculation progresses to refine the fit of the snake. The energy functional using scale space continuation is:

$$E_{edge} = -|G_{\sigma} * \nabla^2 I|^2$$

Where,  $G_{\sigma}$  is a Gaussian with standard deviation  $\sigma$ . Minima of this function fall on the zero-crossings of  $G_{\sigma} \nabla^2 I$ , which define edges as per Marr-Hildreth theory [3].

## 6.6 TERMINATION FUNCTIONAL

Curvature of level lines in a slightly smoothed image can be used to detect corners and terminations in an image. Using this method, let  $C(x, y)$  be the image smoothed by:

$$C(x, y) = G_{\sigma} * I(x, y)$$

With gradient angle

$$\theta = \arctan\left(\frac{C_y}{C_x}\right)$$

Unit vectors along the gradient direction

$$\mathbf{n} = (\cos \theta, \sin \theta)$$

And unit vectors perpendicular to the gradient direction

$$\mathbf{n}_\perp = (-\sin \theta, \cos \theta)$$

The termination functional of energy can be represented as:

$$E_{term} = \frac{\partial \theta}{\partial n_\perp} = \frac{\partial^2 C / \partial n_\perp^2}{\partial C / \partial n} = \frac{C_{yy}C_x^2 - 2C_{xy}C_xC_y + C_{xx}C_y^2}{(C_x^2 + C_y^2)^{3/2}}$$

## 6.7 CONSTRAINT ENERGY

Some systems, including the original snakes implementation, allowed for user interaction to guide the snakes, not only in initial placement but also in their energy terms. Such constraint energy  $E_{con}$  can be used to interactively guide the snakes towards or away from particular features.

## 6.8 OPTIMIZATION THROUGH GRADIENT DESCENT

Given an initial guess for a snake, the energy function of the snake is iteratively minimized. Gradient descent minimization is one of the simplest optimizations, which can be used to minimize snake energy. The algorithm takes one step at each iteration in the negative gradient of the point with controlled step size  $\gamma$  to find local minima. This gradient-descent minimization can be implemented as

$$\bar{v}_i \leftarrow \bar{v}_i + F_{snake}(\bar{v}_i)$$

Where,  $F_{snake}(\bar{v}_i)$  is the force on the snake, which is defined by the negative of the gradient of the energy field.

$$F_{snake}(\bar{v}_i) = -\nabla E_{snake}(\bar{v}_i) = -\left(w_{internal}\nabla E_{internal}(\bar{v}_i) + w_{external}\nabla E_{external}(\bar{v}_i)\right)$$

Assuming the weights  $\alpha(s)$  and  $\beta(s)$  are constant with respect to  $s$ , this iterative method can be simplified to

$$\bar{v}_i \leftarrow \bar{v}_i - \gamma \left\{ w_{internal} \left[ \alpha \frac{\partial^2 \bar{v}}{\partial s^2}(\bar{v}_i) + \beta \frac{\partial^4 \bar{v}}{\partial s^4}(\bar{v}_i) \right] + \nabla E_{ext}(\bar{v}_i) \right\}$$

## 6.9 DISCRETE APPROXIMATION

In practice, images have finite resolution and can only be integrated over finite time steps  $\tau$ . As such, discrete approximations must be made for practical implementation of snakes. The energy function of the snake can be approximated using the discrete points on the snake.

$$E_{snake}^* \approx \sum_1^n E_{snake}(\bar{v}_i)$$

Consequently, the forces of the snake can be approximated as

$$F_{snake}^* \approx -\sum_1^n \nabla E_{snake}(\bar{v}_i)$$

Gradient approximation can be done through any finite approximation method with respect to  $s$ , such as Finite difference.

## 7. PROCEDURE

The GUI of the processing application was developed in MATLAB R2017a revision on HP PC in Microsoft windows 10 pro OS, with Intel(R) Processor Core(TM) i7-8565U CPU @ 1.80GHz, 1992 MHz, 4 Core(s), and 8 Logical processor(s). Standard MATLAB's GUIDE design tool was used to design the appearance of the GUI. To run the application, either you will run **MIA\_App.fig** in the MATLAB command line or use guide tool to load the GUI. Any way, **MIA\_App.fig** is the main GUI file. The `guide` command should be typed in MATLAB command line, then, **MIA\_App.fig** file must be selected from the path where the file is located. Then the current folder root must be changed to the project root (right click on the GUI file and from the top-up menu select "change current folder as root"). The input data are composed of 64 DICOM image files in the Database, which was included in the application root. So if you are planning to run the GUI with different DICOM image, the image folder shall be copied into the root directory of the project. The application has three main panels including:

- **Actions**
- **Snake Contour**
- **3D Representation**

Also, there is a **Status** panel to show the status and number of loaded images in every action. Below, the above panels and their functionality will be described.

### 7.1 ACTIONS

- **Initial status**

The "action" panel includes Load DICOM, Anonymize, Converting from DICOM to JPG, and vice versa, and finally smoothing function as shown in Figure 3. At the beginning, **Anonymize**, **DICOM to JPG** and **Smooth** buttons are disabled, and **Next** and **Previous** buttons are invisible and Load DICOM and JPG to DICOM buttons are active. When Load DICOM function loaded the DICOM images then they will be active as well to be used during the further stages of the process.



Figure 3: actions panel: Load DICOM and JPG to DICOM buttons are active

- **Load DICOM**

Multiple DICOM files are loaded by pushing "Load DICOM" button; then images will be loaded and displayed in the image windows allocated at the right side of the menu of the GUI and can be navigated through **Next** and **Previous** buttons. As soon as user loads DICOM files, the **Next** button will be visible and when navigates to the next image, **Previous** button will be visible. Similarly, for the last image, **Next** button will be invisible and only **Previous** button will be visible. (As explained before, the target Database has already been located in the root directory of the project and is defined as default path for **Load DICOM** function)

The name of each image is shown at the top of the image when it is shown and the **Patient's Information** is displayed at the bottom of the image. Also, the status of process and number of loaded images are displayed in the **Status** panel.

- **Anonymize**

After selecting DICOM files, the **Anonymize** button is enabled, and by clicking on this button, PatientName, PatientID, and BirthDate fields are anonymized; images can be saved in the selected path and the completion of the operation will be notified using a message as shown in Figure 4.

**Success: Annonimized images have been saved.**

Figure 4: successful anonymized action message

If the images have not been anonymized correctly it can be observed when loading these files.

- **DICOM to JPG**

User can convert DICOM images into JPG format just by clicking on the **DICOM to JPG** button and save the JPG file and extracted **mat** file in the selected path. The result of this action is \*.JPG and \*.mat files which includes DICOM metadata. Finally, the below message will be displayed to inform the user that the conversion was done successfully and the resulted images have been successfully stored in the root directly and ready for further processes.

**Success: DICOM images have been saved as JPG.**

Figure 5: successful DICOM to JPG message

- **JPG to DICOM**

“JPG to DICOM” function performs the reverse conversion of the previous step. The JPG image and related metadata must be selected to allow the application to convert JPG to DICOM. Otherwise, an error message will be shown in Figure 6 to warn the user that one data or both data are missing.

**Error: Please select both jpg and related mat file.**

Figure 6: error message in converting JPG to DICOM

If conversion is done successfully, the completion message will be shown to inform the user as shown in Figure 7.

**Success: JPG images have been saved as DICOM.**

Figure 7: successful JPG to DICOM message

- **Smooth**

A Gaussian filter will be applied to the image selecting a sigma parameter as input value. User should enter the sigma ( $\sigma$ ) value in the dialog box that the default value is 1 as is shown in Figure 8. Smoothing the image is done before applying Snake.




Figure 8: sigma value for Gaussian filter

## 7.2 SNAKE CONTOUR

Snake Segmentation is developed for contour detection as detailed in section 6. After loading DICOM files, **Region** popup menu will be enabled which contains Central Zone (CZ), Transitional Zone (TZ), Peripheral Zone (PZ), Tumour and can be selected by clicking on “--Select Zone—” as is visualized in Figure 9. (After performing segmentation or clicking on the **Next&Previous** button, this component will be reset to the default value.)

The image shows a software panel titled "Snake Contour". It contains a dropdown menu for "Region" with the text "--Select Zone--". Below this are several input fields for parameters: Alpha (0.40), Beta (0.20), Gamma (1.00), Kappa (0.15), W(Eline) (0.30), W(Edge) (0.40), W(Eterm) (0.70), and Iterations (80). At the bottom of the panel are two buttons: "Save" and "Show Segmented Regions".

Figure 9: snake panel

Although default values for snake parameters were defined and hard coded, they are editable. So, just replace the old value with a desire value and click on enter to be replaced by the old value. The definition of each value is as follows (for better understanding of the each value refer section 6)

- $\alpha$  (**alpha**): Elasticity of the snake is defined by this parameter. It controls the tension of the contour by combining with the first derivative
- $\beta$  (**beta**): Rigidity in contour is defined with this, through combination with the second derivative
- $\gamma$  (**gamma**): This parameter defines the step size
- $k$  (**kappa**): A scaling factor for energy
- $W(E_{line})$ : A weighing factor for intensity based potential
- $W(E_{edge})$ : Weighing factor for edge based potential
- $W(E_{term})$ : Weighing factor for termination potential

The user can define the number of iteration as well by changing the default value of "iteration" parameter.

To perform the snake segmentation, first, one of the regions is selected just by clicking on the **Region** and selecting one option from popup menu. Selecting the initial position of the snake contour is done through clicking on the image and selecting control points, which later, are interpolated (Spline based) into a contour. Also the selected images can be navigated through **Next** and **Previous** buttons to find an image, which has all of the zones.

Before segmentation, **Save** and **Show Segmented Regions** are disabled. After segmentation of each zone, 3 operations are done as shown in Figure 10:

1. Surface and Volume calculation will be done automatically and the results will be displayed in the **Measures** panel.
2. **Save** button will be enabled and the result can be saved as mask in the user-selected path. These masks are used in the **3D Representation** part. (There is a predefined **Segmentation** folder in the project root, which is the default path for saving masks. Depend on the selected zone, suggested name to save is concatenation of the abbreviation zone and selected image name)
3. **Show Segmented Regions** button will be enabled too and will show all of the segmented results together.



### Snake Contour

Region

Central Zone (CZ)

Alpha

0.40

W(Eline)

0.30

Beta

0.20

W(Edge)

0.40

Gamma

1.00

W(Etern)

0.70

Kappa

0.15

Iterations

80

Save

Show Segmented Regions

### Measures

	Surface	Volume
CZ	721.2031	901.5039
TZ		
PZ		
Tumor		

Figure 10: snake contour

### 7.3 3D REPRESENTATION

The panel shown in Figure 11 is responsible for 3D reconstruction, which contains two buttons: **Load Segmented Images** and **Show 3D Representation**. Clicking on the **load segmented images** enables the **show 3D reconstruction** button. The segmented masks for each zone, which were created during the **Snake Contour** processing can be selected by the user using the **load segmented images** and reconstruct 3D.

### 3D Representation

Load Segmented Images

Show 3-D Representation

Figure 11: 3D representation panel

If the selected directory is out of JPG file, the bellow message will be shown.

**Warning: No JPG file has been found.**

Figure 12: warning folder out of JPG file

Just clicking on the **“Show 3D Representation”** button will show the 3D representation of the masks.

## 8. FLOW CHARTS AND ALGORITHMS

The actions were described in section 7 shows only the functionalities that were request in sub-section 2.1. Now, the question is “to what base these functions were designed?” this section will describe how the function were designed and have been implemented. In the next section, one DICOM file from the attached Database (see the attached zip file) is selected and feed in each function as an input data the actions results are generated.

### 8.1 FUNCTION FLOW CHARTS

In Loading DICOM action, the files are loaded and presented in the right figure and patient’s information will be shown in the bottom panel as detailed in section 7. The `dicomread` and `dicominfo` MATLAB methods are used to implement this functionality. Figure 13 shows the flowchart of this action.




Figure 13: load DICOM function's flow

The **dicomanon** MATLAB function is used to anonymize the selected DICOM files. The related flow is illustrated in Figure 14.




Figure 14: Anonymize function's flow

The DICOM metadata will be saved separately as **.mat** file in the converting DICOM to JPG action. In Figure 15 the related flow is shown.




Figure 15: DICOM to JPG function's flow

In the opposite action, i.e. JPG to DICOM, JPG file and related .mat file will be merged together.




Figure 16: JPG to DICOM function's flow

Applying Gaussian filter was developed using **filter2** and **fspecial** function using  $\sigma$  as paramet Figure 17Error!  
Reference source not found.Error! Reference source not found..




Figure 17: smooth function's flow

Segmentation on the objective zones including Central Zone, Transitional, Peripheral and Tumour have been done one by one according to flow in Figure 18 and using theory described in section 6 and parameter setting as detail in section 7.




Figure 18: Snake function's flow

After selecting the ROI points, system will automatically start making our snake move based on the ROI selected and parameters given. This has been achieved using below function:

```
iterate(image, xs, ys, alpha, beta, gamma, kappa, wl, we, wt, iterations, color)
```

Two functions `calculateSurface` and `calculateVolume` are developed in this part. To calculate the surface, using `regionprops` MATLAB function, we can get the pixel area and sum of them multiply `PixelSpacing` gets final surface. Also volume equals to surface multiply by `sliceThickness`. Figure 19 shows both of their flows.




Figure 19: (a) calculate Surface's function flow (b) calculate volume's function flow

The 3D reconstruction uses the segmented images and puts them into an arranged order as of original order. Then uses surface function to plot them. Then the display parameters are set in another function as shown in Figure 20.




Figure 20: 3D representation function's flow

## 9. RESULTS AND FUNCTIONAL DISPLAY

In this section, result of the above functions will be presented and detail. All of the images of a given patient are selected to load as is shown in Figure 21. As seen, the image is loaded in the image windows, the image file name in the directly is displayed at the top and the patient's information is displayed at the bottom of the image. There are already some images before this image because next and previous buttons are active to allow the user to go backward and forward to see more photos. The results of the anonymized DICOM files will be shown that the objective fields, i.e. PatientName, PatientID, and BirthDate are modified.




Figure 21: DICOM files

The status and number of loaded images are displayed in Status panel (Figure 22Figure 22Figure 22).




Figure 22: status panel




Figure 23: anonymized DICOM file

The image is converted from DICOM to JPG and MAT format and stored in the root directory of the project. So, as seen in Figure 24, for instance for given patient, image00002.jpg and image00002.mat have been generated after clicking on DICOM to JPG button.






Figure 24: JPG and mat files

Gaussian filtering with  $\sigma = 1$ , on a given image has produced the result as shown in Figure 25.




Figure 25: smoothed image

The results are illustrated in Figure 26 and related snake parameters are shown in Figure 27. Surface area and volume of each zone will be shown in Figure 28.




Figure 26: segmented DICOM file (a) peripheral zone (b) transitional zone (c) tumour (d) central zone

Snake Contour

Region

Tumor

Alpha

0.40

W(Eline)

0.30

Beta

0.20

W(Edge)

0.40

Gamma

1.00

W(Eterm)

0.70

Kappa

0.15

Iterations

80

Figure 27: snake parameters

By pushing the save button, **mask** and **snake** result, as the same name of selected image, is saved. The default path has been set to **Segmentation/masks** and **Segmentation/snake** folders in the root. The selected zones the masks are created by segmentation are visualized in Figure 29Figure 29.

Measures

	Surface	Volume
CZ	972.9141	1216.1426
TZ	1288.5391	1610.6738
PZ	687.2734	859.0918
Tumor	278.5391	348.1738

Figure 28: surface and volume of each zone



Figure 29: created masks

Figure 30 shows all segmented zones in one image.




Figure 30: all zones in one image

As mentioned in Section 8.1, these masks are used as input data for 3D representation. 3D representation of zones and full gland respectively are shown in Figure 31 and Figure 32.




Figure 31: 3D representation of zones




Figure 32: 3D representation of full prostate gland

## 10. CONCLUSION

A tool for prostate DICOM image analysis were presented and different functionalities were detailed. Among these functionalities, loading DICOM images and converting them into JPG and vice versa as well as some preprocessing, notably smoothing, were presented. The snake segmentation was detail and the volume and surface extraction as well. The flow chart behind each function was presented and the implementation was explained. At the end some result were presented to show how the tool could be employed for prostate image

analysis. This tool is a useful tool for doctor to load and manipulate DICOM image and study each patient for cancerous cell diagnosis.

## 11. REFERENCES

- [1] Anil Bhavsar and SadhnaVerma. Anatomic imaging of the prostate. *BioMed research international*, 2014, 2014.
- [2] Witkin A. Terzopoulos D. Kass, M. Snakes: active contourmodels. *Int. J. Comput. Vis.*1(14), 321â€“331, 1998.
- [3] David Marr and Ellen Hildreth. Theory of edge detection. *Proceedings of the Royal Society of London. Series B. Biological Sciences*, 207(1167):187–217, 1980.
- [4] Gillian Murphy, MasoomHaider, SangeetGhai, and BoraiahSreeharsha. The expanding role of mri in prostate cancer. *American Journal of Roentgenology*, 201(6):1229–1238, 2013.
- [5] Turkbey, Baris, and Peter L. Choyke. "Multiparametric MRI and prostate cancer diagnosis and risk stratification." *Current opinion in urology* 22.4 (2012): 310.
- [6] Hegde, John V., et al. "Multi-parametric MRI of prostate cancer: An update on state-of-the-art techniques and their performance in detecting and localizing prostate cancer." *Journal of Magnetic Resonance Imaging* 37.5 (2013): 1035-1054.

## 12. APPENDIX

### 12.1 LOADDICOM

```
selectedFile = getSelectedDICOMFile('./Database/*','MultiSelect')
patientInformation = dicominfo(selectedFile)
image = dicomread(patientInformation)
show(image)
show(patientInformation )
```

### 12.2 ANONYMIZE

```
writePath = selectWritePath();
fori=1 to loadedImageCounts:
oldFile = concat(path,images(i))
newFile = concat(writePath ,images(i))
dicomanon(oldFile,newFile,'keep',
('StudyID','StudyDate','SliceLocation','InstanceNumber'))
end
```

### 12.3 DICOMTOJPG

```
writePath = selectWritePath();
fori=1 to loadedImageCounts:
D = dicominfo(file_name);
%extract the metadata from dicom image
matfile = fullfile(write_path,concat(filename,'.mat'));
%save each image's meta data as .mat file
save (matfile);
%extract the image from dicom
dcmImage = dicomread(D);
%save it to jpeg file
dcmImagei = convertToUint8Format(dcmImage);
imwrite(dcmImagei, concat(write_path,file_name,'.jpg'), 'jpg');
end
```

## 12.4 JPGTODICOM

```
names = getSelectedJPGMatfiles();
allMeta = dir(matPath);
allImage = dir(imagePath);
for i=1 to length(names)/2:
    %let's get the JPEG image
    I = imread(completeImageName);
    currentMeta = load(completeMatName);
    %write the DICOM image in desired directory with metadata
    dicomwrite(I, currentMeta.D);
end
```

## 12.5 SMOOTH

```
sigma = getSigma();
If sigma > 0
    smask = fspecial('gaussian', ceil(3*sigma), sigma);
    image = filter2(smask, image, 'same');
    showImage();
```

## 12.6 SNAKE

```
%implementation of "Snakes: Active Contour Models"
[xs, ys] = getSnake(image)
% Let's make the snake move
[region_xs, region_ys] = interate(image, xs, ys, alpha, beta_val, gamma,
kappa, weline, weedge, weterm, iteration, 'r');
%create mask
Io=roipoly(image,xs,ys);
%find the surface & volume
pixel_area = calculateSurface(Io,info.PixelSpacing);
pixel_volume = calculateVolume(pixel_area,info.SliceThickness);
```

## 12.7 CALCULATESURFACE

```
%this function is used to calculate area
surface = calculateSurface(image,pixelSpacing)
stats = regionprops(image, 'Area'); %to get region properties
allSurfaces = [stats.Area]; %to get all areas
surface = sum(allSurfaces); % sum up all pixels
surface = surface * pixelSpacing; % pixels * spacing
```

## 12.8 CALCULATEVOLUME

```
%this function is used to calculate volume using area and slicethickness
volumeImages = calculateVolume( area, sliceThickness );
volumeImages = area * sliceThickness;
```

## 12.9 3D REPRESENTATION

```
%smooth the slices using gaussian filter
J = smooth3(images,'gaussian',5);
% formisosurface and extract patch
p = patch(isosurface(images,100));
% to find isonormals
isonormals(images,p);
%set the properties of patch and draw it
set(p, 'FaceColor', 'red', 'EdgeColor', 'none', 'FaceAlpha', '0.6');
draw();
```